

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS

TITAN  
FILTRATION SYSTEMS



SUN SRI



TITAN

FILTRATION SYSTEMS



## Dissolution and HPLC Syringe Filters

### CONSISTENT RESULTS, CONSISTENT VALUE

Filtration is a critical and valuable strategy to ensure consistent and reliable experimental procedures. From method development investigation to validated quality control techniques, the superior quality of **TITAN** Filtration Products from *sun sri*, secures your results.

### Titan2 Syringe Filters


- Standard 1µm Prefilter in Our 30mm Filters Ensure Increased Sample Handling Capacity
- Sample Distribution Rings Maximize Filtration Area Uniformity and Speed
- Low Extractable Levels Generate Interference Free Data
- Unequivocal Reproducibility, Filter to Filter, Lot to Lot.

### DESIGNED FOR RESULTS, ENGINEERED FOR PERFORMANCE

Not just concerned with results, *Titan2* Filters are designed and produced with exceptional features that assure un-interrupted performance. *sun sri* has addressed the most common problems experienced with syringe filters and not only solved these issues, but guarantees *Titan2* performance, without exception.

- Pressure Rated to 90 psi for Unsurpassed Integrity
- Sealing Ring Minimizes Bursting and Leakage
- Chemically Resistant Polypropylene Housings

## TITAN2 ...The Logical Step in Obtaining Superior Results.



*Titan2's color-code ID System permits easy identification of membrane type and porosity. Our unique sealing ring is coded for membrane type (i.e., Nylon, PTFE, PVDF, etc.) and porosity (i.e., 0.45µm, 0.2µm, 5.0µm, etc.) for instant visual verification. The Titan System goes beyond the filter itself and includes our re-usable, rigid, transparent, color-coded containers. Now maintain a proper inventory of the filter you need at a moments glance.*



## Titan2 Dissolution Testing Syringe Filters

Regulatory requirements demand that validated methods be robust and rugged to withstand operator dependent procedures and product variability. Titan Filters assure consistency and reproducibility with the strict quality assurance production methodology of sun sri.



## Titan2 FILTERS ARE QUALITY MONITORED FOR:

- Low UV Extractable for Superior Results in Spectrophotometric or HPLC Procedures
- All Critical Filter Specifications are Tested and Confirmed... Bubble Point, Burst Pressure, Membrane Retention/Exclusion, Flow Rate, UV Extractables
- Retention Samples Kept for Six Months on All Lot Production Permits sun sri to Validate and Assist in Solving YOUR Problem
- Actual HPLC Performance and Reproducibility Validation by a Standard Test Mixture Run.

	Transp. Blue PTFE 0.2µm
	Light Brown Reg. Cellulose 0.45µm
	Dark Brown PTFE 5.0µm
	Amber Glass Microfibre 0.7µm
	Navy Blue PTFE 1.0µm
	Royal Blue Polypropylene 0.45µm
	Orange-Yellow Polysulfone 0.45µm
	Opaque Med. Yellow Nylon 1.2µm
	Transp. Yellow PTFE 0.45µm
	Transp. Green Nylon 0.45µm
	Light Green Cellulose Acetate 0.20µm
	Pink Nylon 5.0µm
	Light Purple Nylon 0.20µm
	Dark Purple Glass Microfibre 3.1µm
	Transp. Orange Cellulose Acetate 0.45µm
	Orange Glass Microfibre 1.2µm
	BLACK PVDF 0.2µm
	Dk. Gray Polysulfone 0.20µm
	Red PVDF 0.45µm
	White Nylon 0.45µm
	Granite Reg. Cellulose 0.20µm

## Nylon Syringe Filters

### Membrane of Choice for Analytical Applications

- Naturally hydrophilic with broad chemical resistance
- Excellent flow rates and high throughput loading
- Standard 1µm prefilter (30mm) for optimal performance with particulate-laden analytes
- Certified for chromatographic performance with an actual HPLC analytical standards run

Titan2 Nylon Syringe Filters offer universal application for analytical procedures. Hydrophilic Nylon is extremely well suited for aqueous or organic sample prep and HPLC, GC or dissolution sample analysis. Due to its excellent flow characteristics and mechanical stability, Nylon offers the best combination of physical parameters to meet the most stringent analytical needs.



#### Specifications:

Membrane: HPLC Certified Nylon	Max. Operating Temp.	100°C
Housing: Medical Grade, Virgin Polypropylene	Max. Operating Pressure:	4mm - 75 psi
Prefilter: Binder Free Glass Microfiber		17mm - 115 psi
		30mm - 90 psi
Connections: Enhanced Female Luer-Lok Inlet	Retention Volumes:	4mm - <15 µL
Male Slip Outlet		17mm - <29 µL
		30mm - <137 µL

Autoclave: Sterilize by dry heat at 121°C for 15 min.

#### Applications:

HPLC and GC Sample Prep and clean up  
Dissolution Sample Analysis  
General Sample Prep prior to Analytical Analysis  
Mixed sample matrix of aqueous or organic dissolved analytes.


#### Chemical Incompatibilities:

Acids >1N  
Halogenated solvents  
Proteinaceous Samples with high non-specific binding affinities

#### ORDERING INFORMATION

Identifying Color Code	Porosity	Std 1µm Prefilter	Catalog Number	Recommended Sample Vol.	Qty
<b>4mm</b>					
LT PURPLE	0.20µm	N	42204-NN	<2mL	100/pkg
TRANSF GREEN	0.45µm	N	44504-NN	<2mL	100/pkg
<b>17mm</b>					
LT PURPLE	0.20µm	N	42213-NN	<30mL	200/pkg
TRANSF GREEN	0.45µm	N	44513-NN	<30mL	200/pkg
<b>30mm</b>					
LT PURPLE	0.20µm	Y	42225-MN	<150mL	100/pkg
	0.20µm	Y	42225-NC	<150mL	500/cs
TRANSF GREEN	0.45µm	Y	44525-MN	<150mL	100/pkg
	0.45µm	Y	44525-NC	<150mL	500/cs
WHITE	0.45µm	N	44526-MN	<100mL	100/pkg
TRANSF GREEN	1.2µm	N	41225-MN	<100mL	100/pkg
PINK	5.0µm	N	45025-MN	<100mL	100/pkg

#### Throughput Volume Comparison- Titan & Titan2 Nylon Filters



Titan2 offers a considerable enhancement to throughput volume with our new filter improvements compared to our first generation Titan syringe filters. Although our original, pre-filtered Titan Syringe Filters were already an optimal choice over competitive filters, our new Titan2 Syringe Filters bring throughput value differences to an even greater level.

## PTFE Syringe Filters

*Excellent chemical resistance for use with organic matrices*

- Standard 1µm prefilter (30mm) for optimal performance with particulate-laden samples
- Naturally hydrophobic membrane provides excellent flow rates and high loading
- Certified for chromatographic performance with an actual HPLC analytical standards run
- Exceptional temperature stability provides extended sampling range over other membranes



Titan2 PTFE Syringe Filters are indicated for filtration of gaseous or organic solvent based samples. The natural hydrophobic membrane exhibits broad chemical resistance and unsurpassed temperature stability to address aggressive sample matrixes and extreme temperature situations. Due to the membranes hydrophobic nature, PTFE filters can be utilized as a moisture barrier in venting applications.

### Specifications:

Membrane: HPLC Certified PTFE,  
w/polypropylene support

Housing: Medical Grade, Virgin Polypropylene

Prefilter: Binder Free Glass Microfiber Prefilter

Connections: Enhanced Female Luer-Lok Inlet  
Male Luer Slip Outlet

Max. Operating Temp. 100°C

Max. Operating Pressure: 4mm - 75 psi

17mm - 115 psi

30mm - 90 psi

Retention Volumes: 4mm - <15 µL

17mm - <29 µL

30mm - <137 µL

Autoclave: Sterilize by dry heat at 121°C for 15 min.

### Applications:

HPLC and GC Sample Prep and clean up  
Dissolution Sample Analysis-HCl media based  
General Sample Prep prior to Analytical  
Analysis

Elevated temperature samples, caustic or  
acidic solutions


### Chemical Incompatibilities:

Perchloric Acid  
Methylene Chloride (limited exposure)  
Dioxane, DMF, Formic Acid >50%  
Aqueous based sample matrix (unless filter is  
pre-wetted with an alcohol)

### ORDERING INFORMATION

Identifying Color Code	Porosity	Std 1µm Prefilter	Catalog Number	Recommended Sample Vol.	Qty
<b>4mm</b>					
ROYAL BLUE	0.20µm	N	42204-NP	<2mL	100/pkg
TRANSF YELLOW	0.45µm	N	44504-NP	<2mL	100/pkg
<b>17mm</b>					
ROYAL BLUE	0.20µm	N	42213-NP	<30mL	200/pkg
TRANSF YELLOW	0.45µm	N	44513-NP	<30mL	200/pkg
<b>30mm</b>					
ROYAL BLUE	0.20µm	Y	42225-NMP	<150mL	100/pkg
ROYAL BLUE	0.20µm	Y	42225-NPC	<150mL	500/cs
TRANSF YELLOW	0.45µm	Y	44525-NP	<150mL	100/pkg
TRANSF YELLOW	0.45µm	Y	44525-PC	<150mL	500/cs
ROYAL BLUE	1.0µm	N	41025-NP	<100mL	100/pkg
OK BROWN	5.0µm	N	45025-NP	<100mL	100/pkg

### Validated HPLC Performance- Caffeine (µg/mL)



Only Titan2 Filters are validated for HPLC performance and reproducibility by an actual HPLC standards run. Each lot is quality controlled for retention time, peak area, and peak shape of our 100% inorganic caffeine standard. This assures you of outstanding performance, reproducibility and zero eluate/leak conditions. For other filter specifications order actual run conditions for Titan2.

## PVDF Syringe Filters

*Low non-specific binding with exceptional chemical resistance*

- Standard 1µm prefilter (30mm) for optimal performance with particulate-laden samples
- Hydrophilic membrane provides excellent flow rates and low binding coefficients
- Certified for chromatographic performance with an actual HPLC analytical standards run
- Compatible with a wide range of aqueous and organic based sample environments



### Specifications:

Membrane: HPLC Certified PVDF  
Housing: Medical Grade, Virgin Polypropylene  
Prefilter: Binder-Free Glass Microfiber Prefilter

Connections: Enhanced Female Luer-Lok Inlet  
Male Luer Slip Outlet

Autoclave: Sterilize by dry heat at 121°C for 15 min.

Max. Operating Temp.: 100°C  
Max. Operating Pressure: 4mm - 75 psi  
17mm - 115 psi  
30mm - 90 psi  
Retention Volumes: 4mm - <15 µL  
17mm - <29 µL  
30mm - <137 µL

### Applications:

HPLC and GC Sample Prep and clean up  
Protein based samples with high non-specific binding

### Chemical Incompatibilities:

DMF, DMSO, MEK, Acetone and most caustic solutions >6N

### ORDERING INFORMATION

Identifying Color Code	Porosity	Std 1µm Prefilter	Catalog Number	Recommended Sample Vol.	Qty
<b>4mm</b>					
BLACK	0.20µm	N	42204-PV	<2mL	100/pkg
RED	0.45µm	N	44501-PV	<2mL	100/pkg
<b>17mm</b>					
BLACK	0.20µm	Y	42213-PV	<30mL	200/pkg
RED	0.45µm	N	44513-PV	<30mL	200/pkg
<b>30mm</b>					
BLACK	0.20µm	Y	42225-PV	<150mL	100/pkg
RED	0.45µm	N	44525-PV	<150mL	100/pkg

*Please go to page 11 for full information about our All-Plastic Disposable Syringes*



## Glass Microfiber (GMF) Syringe Filters

*Low extractable, 100% binder-free, borosilicate glass filter for dissolution analysis*

- Increased sample throughput with 1.2 µm and 3.1 µm porosity filters
- Excellent for use with viscous or particle laden samples coupled with broad range chemical compatibility

### Specifications:

Membrane: Binder-free glass microfiber  
Housing: Medical Grade, Virgin Polypropylene  
Connections: Enhanced Female Luer-Lok Inlet  
Male Luer Slip Outlet

Max. Operating Temp.: 110°C  
Max. Operating Pressure: 30mm - 90 psi  
Retention Volumes: 30mm - <137 µL

### Applications:

Viscous or particle laden samples,  
Dissolution testing

### Chemical Incompatibilities:

Limited resistance with Ammonium, Potassium, and Sodium Hydroxide Solutions

### ORDERING INFORMATION

Identifying Color Code	Porosity	Catalog Number	Recommended Sample Vol.	Qty
<b>30mm</b>				
BROWN	0.70µm	40725-GM	<150mL	100/pkg
ORANGE	1.20µm	41225-GM	<175mL	100/pkg
DE PURPLE	3.10µm	42725-GM	<175mL	100/pkg

## Regenerated Cellulose Syringe Filters


*Superior chemical resistance optimized for biological sample recoveries*

- Hydrophilic membrane provides excellent flow rates and extremely low binding coefficients
- Superior choice for biological assays, gel capsule dissolution testing, protein sample matrices
- Compatible with a wide range of aqueous and organic based sample environments
- Certified for chromatographic performance with an actual HPLC analytical standards run



Titan2 Regenerated Cellulose Syringe Filters are ideally suited for almost any laboratory procedure from HPLC sample prep to dissolution sample validation. Regenerated Cellulose possesses superior chemical resistance in either aqueous or organic sample environments. Its extremely low biological based binding coefficient is ideally suited for maximum sample recoveries of biological based assays. Regenerated Cellulose contains no binders, surfactants or wetting agents to assure zero extractables in analytical procedures.

BINDING CHARACTERISTICS of VARIOUS MEMBRANES  
(Binding in  $\mu\text{g/ml}$  of G-Gobulin)



### Specifications:

Membrane: HPLC Certified Regenerated Cellulose	Max. Operating Temp.	100°C
Housing: Medical Grade, Virgin Polypropylene	Max. Operating Pressure:	4mm - 75 psi 17mm - 115 psi 30mm - 90 psi
Protein Binding: <math><5\mu\text{g}/\text{cm}^2</math>	Retention Volumes:	4mm - <math><15\ \mu\text{L}</math> 17mm - <math><29\ \mu\text{L}</math> 30mm - <math><137\ \mu\text{L}</math>
Connections: Enhanced Female Luer-Lok Inlet Male Luer Slip Outlet		
Autoclave: Sterilize by dry heat at 121°C for 15 min		

### Applications:

HPLC and GC Sample Prep and clean up  
Dissolution sample assays, especially high binding tablets or capsules  
Protein based samples with high non-specific binding  
Sample analysis with required maximum recoveries  
Analysis requiring low non specific binding over a wide pH range

### Chemical Incompatibilities:

Sulfuric Acid, Hydrochloric Acid, Phosphoric Acid or Nitric Acid >25%, DMF, Phenol

### ORDERING INFORMATION

Identifying Color Code	Porosity	Std 1um Prefilter	Catalog Number	Recommended Sample Vol.	Qty
<b>4mm</b>					
WHITE	0.20 $\mu\text{m}$	N	52204-RC	<math><2\text{mL}</math>	100/pkg
LIGHT BROWN	0.45 $\mu\text{m}$	N	54504-RC	<math><2\text{mL}</math>	100/pkg
<b>17mm</b>					
WHITE	0.20 $\mu\text{m}$	Y	52213-RC	<math><30\text{mL}</math>	200/pkg
LIGHT BROWN	0.45 $\mu\text{m}$	N	54513-RC	<math><30\text{mL}</math>	200/pkg
<b>30mm</b>					
WHITE	0.20 $\mu\text{m}$	Y	52225-RC	<math><150\text{mL}</math>	100/pkg
LIGHT BROWN	0.45 $\mu\text{m}$	N	54525-RC	<math><150\text{mL}</math>	100/pkg

## PES Syringe Filters - Ion Chromatography Certified

Precise results in sensitive analysis of ionic analytes

- Certified for low level IC interference by ICP analysis
- Hydrophilic membrane provides excellent flow rates and low binding coefficients
- Low affinity for binding drugs; ideal dissolution testing procedures



### Specifications:

Membrane: ICP Certified PES (PolyEtherSulfone)  
Housing: Medical Grade, Virgin Polypropylene

Connections: Enhanced Female Luer-Lok Inlet  
Male Slip Outlet

Autoclave: Sterilize by dry heat at 121°C for 15 min.

Max. Operating Temp. 100°C  
Max. Operating Pressure: 4mm - 75 psi  
17mm - 115 psi  
30mm - 90 psi  
Retention Volumes: 4mm - <15 µL  
17mm - <29 µL  
30mm - <137 µL

### Applications:

IC (ion chromatography) sample prep & analysis  
Dissolution testing  
Protein based sample in aqueous solutions

### Chemical Incompatibilities:

Concentrated Acids, MeCl, Chloroform, Hexane, Acetone, MEK, THF, DMSO

### ORDERING INFORMATION

Identifying Color Code	Porosity	Std 1µm Prefilter	Catalog Number	Recommended Sample Vol.	Qty
<b>17mm</b>					
	0.20µm	N	42213-PS	<2mL	200/pkg
	0.45µm	N	44513-PS	<2mL	200/pkg
<b>30mm</b>					
	0.20µm	Y	42225-PS	<30mL	100/pkg
	0.45µm	N	44525-PS	<30mL	100/pkg

Analyte	Method Detection Level (ppb)	Analysis Results (ppb)	
		0.2µm PES	0.45µm PES
Bromide	<100	< MCL	< MCL
Chloride	20	< MCL	< MCL
Fluoride	20	< MCL	< MCL
Nitrate as N	50	< MCL	< MCL
Oxalacetate as P	<100	< MCL	< MCL
Sulfate	<100	< MCL	< MCL
Borates	2	< MCL	< MCL
Calcium	10	30.9	38.9
Phosphate	50	80	70
Magnesium	20	< MCL	< MCL
Sodium	20	80.8	80

## Polypropylene Syringe Filters

Chemically resistive membrane with low protein binding

- Hydrophilic membrane for aqueous or organic sample matrixes
- Use with protein or peptide based assays



### Specifications:

Membrane: Hydrophilic Polypropylene  
Housing: Medical Grade, Virgin Polypropylene

Connections: Enhanced Female Luer-Lok Inlet  
Male Slip Outlet

Max. Operating Temp. 110°C  
Max. Operating Pressure: 17mm - 115 psi  
30mm - 90 psi  
Retention Volumes: 17mm - <29 µL  
30mm - <137 µL

### Applications:

Protein or peptide based assays, general HPLC analysis

### Chemical Incompatibilities:

Hexane, Toluene, Benzene, limited resistance to MeCl and Chloroform

### ORDERING INFORMATION

Identifying Color Code	Porosity	Catalog Number	Recommended Sample Vol.	Qty
<b>17mm</b>				
	0.20µm	42213-PP	<30mL	200/pkg
	0.45µm	44513-PP	<30mL	200/pkg
<b>30mm</b>				
	0.20µm	42225-PP	<150mL	100/pkg
	0.45µm	44525-PP	<150mL	100/pkg

## Cellulose Acetate Syringe Filters

For filtering of aqueous solutions or biological samples

- Hydrophilic membrane provides excellent flow rates and extremely low testing, protein binding coefficients
- Superior choice for biological assays, gel capsule dissolution testing, protein sample matrices
- Certified for chromatographic performance with an actual HPLC analytical standards run



### Specifications:

Membrane: HPLC Certified Cellulose Acetate  
Housing: Medical Grade, Virgin Polypropylene  
Protein Binding: <24µg/cm<sup>2</sup>

Connections: Enhanced Female Luer-Lok Inlet  
Male Luer Slip Outlet

Autoclave: Sterilize by dry heat at 121°C for 15 min.

Max. Operating Temp. 110°C  
Max. Operating Pressure: 4mm - 75 psi  
17mm - 115 psi  
30mm - 90 psi  
Retention Volumes: 4mm - <15 µL  
17mm - <29 µL  
30mm - <137 µL







### Applications:

Protein based samples with high non-specific binding  
Sample analysis with required maximum recoveries

### Chemical Incompatibilities:

Acids, NaOH, Dichloromethane, Chloroform, Ketones, DMSO, THF

### ORDERING INFORMATION

Identifying Color Code	Porosity	Catalog Number	Recommended Sample Vol.	Qty
<b>4mm</b>				
	0.20µm	42204-CA	<2mL	100/pkg
	0.45µm	44504-CA	<2mL	100/pkg
<b>17mm</b>				
	0.20µm	42213-CA	<30mL	200/pkg
	0.45µm	44513-CA	<30mL	200/pkg
<b>30mm</b>				
	0.20µm	42225-CA	<150mL	100/pkg
	0.45µm	44525-CA	<150mL	100/pkg

## Titan-MSF™

Microsample, centrifuge filters for rapid particulate removal

- Filter multiple samples simultaneously
- Low extractable device with maximum sample recoveries
- Compatible with a wide range of aqueous and organic based sample environments

Titan MSF microsample filters permit rapid removal of particulates with minimized sample retention (5µL) for efficient sample prep prior to analytical techniques. Using a microcentrifuge, the Titan MSF can process volumes from 200µL to 750µL. The unit is ideal for viscous, colloidal or high solids samples. Units include sample chamber and receiver tube and will fit rotors that accept 1.5 to 2.0mL tubes.



### Specifications:

Membrane: HPLC Certified Nylon & PTFE  
Housing: Medical Grade, Virgin Polypropylene  
Retention Volume: <5µL

Max. Operating RDF: 10,000 x g  
Dimensions: 44.7mm (fully assembled w/cap)  
Max Operating Temp: 80°C

### Applications:

Small volume sample runs for rapid cleanup  
HPLC sample prep, viscous sample matrices  
Removal of precipitates

### ORDERING INFORMATION

Membrane Type	Porosity	Catalog Number	Recommended Sample Vol.	Qty
Nylon	0.20µm	62215-NN	200µL to 750µL	100/pkg
Nylon	0.45µm	64515-NN	200µL to 750µL	100/pkg
PTFE	0.20µm	62215-NP	200µL to 750µL	100/pkg
PTFE	0.45µm	64515-NP	200µL to 750µL	100/pkg



## Power Filtration System

*Optimized speed, power and efficiency in filtration procedures*

- Applies three times greater force compared to manual filtration methods
- Samples that are difficult to filter manually are now easily processed
- Speeds filtration applications compared to manual method, perfect for processing multiple samples
- Certified for chromatographic performance with an actual HPLC analytical standards run

### Laboratory tested and proven

Hot Shot has been laboratory tested by an independent engineering firm to confirm the force and speed advantage the Hot Shot System produces compared to manual filtration with a syringe. Hot Shot was found to produce the following results:

1. Force applied advantage of 3.009 over manually applied pressure
2.  $P_{HotShot} = 87 \text{ psi}$  vs.  $P_{manual} = 25 \text{ psi}$  (5 mL syringe)
3. Elapsed time to fully empty, >6x for test solutions (10mL cartridge)

### A unique system approach, but easy to operate;

The Hot Shot system is built around the permanent filter processor and utilizes disposable filter cartridges, which come preassembled with a 30mm filter attached to either a 5mL or 10mL sample reservoir. Simple to operate, the user fills the preassembled filter cartridge then applies the piston cap. The reservoir piston cap and filter are disposable to eliminate cross contamination between samples. The filled cartridge is then attached to the processor and the trigger is depressed to filter the sample. The units' built-in safety features release the cartridge before the pressure rating of the filter unit is exceeded. The result is no filter bursting or filter "blow-off". The user and lab work area are protected from contamination by sample leakage from over-pressurized filters. Multiple samples are processed quickly without causing user fatigue.



Cartridge Type and Size	Porosity	Catalog No.	Qty
<b>Nylon</b>			
5mL	0.20µm	52205-NH	50/pkg
5mL	0.45µm	54505-NH	50/pkg
10mL	0.20µm	52210-NH	40/pkg
10mL	0.45µm	54510-NH	40/pkg
<b>PTFE</b>			
5mL	0.20µm	52205-HP	50/pkg
5mL	0.45µm	54505-HP	50/pkg
10mL	0.20µm	52210-HP	40/pkg
10mL	0.45µm	54510-HP	40/pkg
<b>Filter Processors</b>			
For 5mL Cartridges	—	59440-HS	1 Ea
For 10mL Cartridges	—	59441-HS	1 Ea

### Intro-Pac Special

*Includes one package of 0.45µm Nylon and PTFE cartridges plus one Hot Shot Filter Processor*

59440-HSP  
5mL Hot Shot Intro Pac Kit

59441-HSP  
10mL Hot Shot Intro Pac Kit